



## Disposable Nitrile Gloves.

1. Made of excellent nitrile synthetic rubber, it avoids the latex-resulting skin allergies.
2. A better performance in puncture resistance, bacterial penetration resistance and chemical resistance.
3. Long service life, flexible operation, soft texture, comfortable wearing experience.
4. Non-poisonous, harmless and odorless. Made with handpicked formula and advanced technologies, the product is soft, comfortable, nonskid and easy to operate, with good protection and physical performances.



<b>Product name :</b>	Disposable Gloves	<b>Payment :</b>	T/T,L/C,etc.
<b>Material :</b>	Nitrile, Powder Free	<b>Grade :</b>	Industrial,Medicaland Food grade
<b>Size :</b>	XS(6-7cm) , S(7-8cm) , M(8-9cm) , L(9-10)cm , XL(10-11cm), the hand width	<b>Packing :</b>	100pcs/box,10boxes/ctn, size:23.4*12*6cm
<b>Color :</b>	Sky blue	<b>Quaiity Standard :</b>	SGS , CNAS , MA , ilac-MRA
<b>Weight :</b>	About 0.042KG/pcs , 0.46KG/box	<b>Certificate :</b>	FDA,CE,ISO,BRC,NSF,TUV,SGS,etc.
<b>MOQ :</b>	10 boxs	<b>Elongation at break(%) :</b>	≥500
<b>Length :</b>	9inch , About 23CM	<b>Tensile Strength(Mpa) :</b>	≥14

**MOQ:** 1000 boxes/10ctns